

Specifications

		ChargePoint® BIO
Containment Performance	Standard	<10µg/m ³
	EXCEL option	<1µg/m ³
Sizes	100mm (4")	●
	150mm (6")	●
	200mm (8")	●
Steam In Place Rating		2 Bar (29 psi)
Product Contact Materials	316L	●
	C22	●
Seals	FKM	●
	FFKM	●
Connection Interface		Hygienic connection to suit process and container
Operation	Manual	●
	Semi Actuated	●
	Fully Actuated	●