

Specifications

		ChargePoint®	ChargePoint® EXCEL	ChargePoint® DL
Containment Performance		<10mcg/m ³	<1mcg/m ³	1-10 mcg/m ³
Size	50mm (2")	●		
	100mm (4")	●	●	
	150mm (6")	●	●	
	200mm (8")	●	●	●
	250mm (10")			●
	300mm (12")			●
Pressure Rating ¹		Up to 6Bar (87psi) ¹	Up to 6Bar (87psi) ¹	None
Product Contact Materials	316L	●	●	●
	Alloy C22	●	●	●
Seals	EPDM	●	●	●
	FKM	●	●	●
	FFKM	●	●	●
Connection Interface ²	Tri-clamp	●	●	Connection to suit IBC/container
	PN6	●	●	
	PN10	●	●	
	ANSI 150#	●	●	
Operation	Manual	●	●	
	Semi Actuated	●	●	
	Fully Actuated	●	●	●

¹ Pressure ratings vary upon the size and version of valve. Active is pressure rated when used in conjunction with a Pressure Rated Plug, Passive or Sightglass. Passive can only be pressure rated when coupled to a pressure rated Active or Pressure Rated Cover.

² Other connection interfaces are available upon request.